



## 1. Introduction

### 1.1. Purpose

## 2. Overview

### 2.1. System Overview

The M20 EEG Device (Component 1 see Fig. 1) has 4 primary components:

- 1) Headset, which consists of an adapter and sensor strip pair,
- 2) Headset accessories: textile caps, sensors, and ear clips/stickers,
- 3) Bioamplifier (bioamp),
- 4) Collect Software [Visualization software]. Note: Collect runs on a Windows laptop referred to as the Visualization laptop




Figure 1: End-to-end M20 System Diagram for use in Clinical Investigation



## 3. References

### 3.1. Applicable Standards

This document demonstrates compliance to the requirements of the following Standards and Regulations:

Document Number	Title
IEC 60086-4	Battery Standards (Lithium Battery)
IEC 60950-1	Power Supply
UL 796	Standard for Printed-Wiring Boards

### 3.2. Internal Documents

The following documents form part of this document. In the event of conflict, this document supersedes:

Document Number	Title

### 3.3. Abbreviations


Acronym/Term	Definition
EEG	Electroencephalography
Task	Individual stimulus test
Session	One or more contiguous Tasks for a single participant



## 4. Architecture

### 4.1. Hardware Architecture


#### 4.1.1. Bioamplifier to / from Adapter (Headset Unit) w/ Sensor Strip



- 4.1.1.1. The Headset Adaptor allows connection of sensors that can be converted the sensor source to the voltage source. The sensors are 2-wire active electrodes that connect to the Adapter via 2-pole touchless connectors. These readings can then be transferred through the connector to the Bioamplifier. The Adaptor includes 9 channels of EEG plus reference and bias connections.



## 4.1.2. Bioamplifier to / from Visualizer PC (w/ collect software system)




- 4.1.2.1. The (4) ADS1299 8 channel bioamplifier and processes received sensor data from the Adapter and passes on to the MCU. The MCU is connected to the ADS1299's and handles the data conversion from the bioamplifiers, communication with the host PC, and monitoring onboard power supplies as well as consolidates EEG data with received 'Mark' data from the STIM source. 2 USB serve as sources for Bioamplifier power. The STIM USB port allows receiving of Stimulus relevant Mark data at the Bioamplifier. The DATA USB port allows receiving of test initiation and execution data from the Visualizer and sends EEG data received to the Visualizer. Input power from 1 of the 2 USB Ports (STIM and DATA) are regulated 'down' by the DC/DC converters. Data sent to the Visualizer is then visualized by the Collect software application continuously.



## 4.2. Software Architecture

### 4.2.1. Collect Software System

#### Collect



- 4.2.1.1. The Collect Software application controls and configures bioamplifier, reads EEG samples from bioamplifier, visualizes EEG data in time and frequency domains for human interpretation
  - 4.2.1.1.1. Records EEG data and metadata to file
  - 4.2.1.1.2. Accepts user input to configure bioamplifier
- 4.2.1.2. Application Framework
  - 4.2.1.2.1. Application with Node.js backend and Javascript/HTML frontend.
  - 4.2.1.2.2. Backend:
    - 4.2.1.2.2.1. Boxer
      - communicates with bioamplifier over serial port
      - writes EEG recording to file.
    - 4.2.1.2.2.2. Serverer
      - Manages IPC (Interprocess Communication) with the Frontend.
  - 4.2.1.2.3. Frontend:



- 4.2.1.2.3.1. Renders UI for EEG metadata input from user, enable recording, change bioamplifier settings
- 4.2.1.2.3.2. Renderer  
Draws EEG data in time domain or frequency domain.
- 4.2.1.2.3.3. Dataer  
Receives EEG data from backend and processes data for drawing.
- 4.2.1.2.3.4. Socketer  
Manages IPC (Interprocess Communication) with the Backend.

## 4.2.2. Bioamplifier Firmware



- 4.2.2.1. Firmware reads EEG samples from the ADS1299's and formats for sending to Visualizer System laptop. Single threaded with each task maintaining a separate state machine. Additionally, the firmware accepts Command Line Interpreter (CLI) commands.
- 4.2.2.2. Firmware utilizes the following primary modules.
  - MAIN
  - EEG
  - SPI
  - LEDS
  - SERIAL
  - CLI
  - EVENTS
  - VERSION
  - SETTINGS
  - Processor Expert Modules (Processor Expert build System) - hardware access



## 5. System Design

### 5.1. Collect Software System

#### 5.1.1. User Interface Design

##### 5.1.1.1. Preset channel display options (dropdown selection)

Configuration selection for channels used for captured EEG data

- M20 (3 ch) - (default) - display of Fz / Cz / Pz channels
- M20 (4 ch) - display of Fz / Cz / Pz / Oz channels
- 16 Channels - 1-16 channels
- 24 Channels - 1-24 channels
- 32 Channels - 1-32 channels
- FSU Wave 1 - Right Earlobe, Left Upper eye, Left Lower Eye, Left Horizontal Eye, Right Horizontal Eye, Fz, Cz, Pz
- FSU Wave 2 - A2, Fp1, Olat, Oinf, Fz, Cz, Pz

##### 5.1.1.2. Site (dropdown selection) - Site where EEG study is performed.

- M20 - (Default)
- FSU
- MSU
- SDSU
- UF
- UU

##### 5.1.1.3. Tasks (drops down selection)

Stimulus tasks associated with captured EEG data

- Doors (Default)
- Emotional Images
- Eyes Closed
- Flankers
- Hot Button
- SRET
- SSVEP
- Oddball
- Test

##### 5.1.1.4. Participant ID (manual entry or arrow selection)

Unique ID for participant to associate captured EEG data

Numeric only

No Default

##### 5.1.1.5. Session (manual entry)

Session number for the participant

Session number default is '1'



- Numeric only
- 5.1.1.6. Port (auto / dropdown selection)  
Selection of port connected to the Bioamplifier
- 5.1.1.7. Record command Button  
Prompts user for saving file and beginning the recording of the collected EEG data. This is synchronized with the execution of Stimulus data execution (sessions) with the participant.
  - 5.1.1.7.1. Save recording to directory prompt  
File directory selection prompt for selecting folder location
  - 5.1.1.7.2. Record button changes to Stop button
  - 5.1.1.7.3. Time (HH:MM:SS) of Record execution elapsed
  - 5.1.1.7.4. Display of file name and directory folder for saved CSV file displayed in prompt when selecting Stop.
  - 5.1.1.7.5. Displays Recording in Graph display
- 5.1.1.8. Configuration Icon  
Display the administrative settings configurable by user.
  - 5.1.1.8.1. Software Version  
Collect software version
  - 5.1.1.8.2. Firmware Version  
Bioamplifier firmware version
  - 5.1.1.8.3. Hardware Version  
Bioamplifier hardware version
  - 5.1.1.8.4. Filter - Order
  - 5.1.1.8.5. Filter - Low Pass - Hz  
Filter Low Pass Hz amount (40 default) - numeric with arrow option
  - 5.1.1.8.6. Filter - High Pass - Hz  
Filter High Pass Hz amount (.1 default) - numeric
  - 5.1.1.8.7. Channels  
Default to selected Preset channels associated  
E.g. M20 (3 Ch)
    - Fz checked
    - Cz checked
    - Pz checkedUser has option to enter channel labels and check which ones they wish to include - alpha numeric
  - 5.1.1.8.8. New Preset Name  
User has the option to enter a new preset name for the selected / entered channels
  - 5.1.1.8.9. Save As New and Apply Command button





- Saves the new preset and adds to the list of Presets available for selection
- 5.1.1.8.10. Send Test Command button
  - Sends a predetermined test message to the Bioamplifier and displays the results on the wave data graph
- 5.1.1.9. Display Layout
  - 5.1.1.9.1. Display Graph display icon selections
    - 5.1.1.9.1.1. (x axis) % x (y axis) Frequency (hz)
    - 5.1.1.9.1.2. (x axis) uV x (y axis) Time (seconds)
  - 5.1.1.9.2. Left hand side (x-axis) Channels
  - 5.1.1.9.3. Bottom (y-axis) numeric Time or Frequency depending on selection.
  - 5.1.1.9.4. Displays Not Recording - Default display
- 5.1.1.10. Filter switch
  - 5.1.1.10.1. Turning on the switch highlights in dark blue and applies the configured filter to the display.
    - 5.1.1.10.1.1. Applies the configured filters
      - Low Pass Hz
      - High Pass Hz was going on yeah what's the plan
  - 5.1.1.10.2. Default is off
- 5.1.1.11. Auto Baseline w/ Continuous checkbox
  - 5.1.1.11.1. Applies a baseline to the graph for each channel
  - 5.1.1.11.2. Default is on
  - 5.1.1.11.3. Continuous checkbox displays the baseline continuously, if not selected the baseline is redrawn repeatedly.
  - 5.1.1.11.4. Not available for the % x Frequency Hz filter
- 5.1.1.12. uV Peak to Peak Amplitude
  - 5.1.1.12.1. Increasing or decreasing options for increasing peak to peak amplitude displayed.
  - 5.1.1.12.2. 1000 uV default
  - 5.1.1.12.3. Not available for the % x Frequency Hz filter
- 5.1.1.13. Waveform data
  - 5.1.1.13.1. Current, continuously received EEG waveform is displayed on the interface from left to right with data display varying depending on configuration and filter settings.
  - 5.1.1.13.2. Mark data received from the stimulus system will be displayed as received (at predetermined intervals that represent stimulus participant actions) along with the associated EEG data.



- 5.1.1.13.3. Each channel configured will display the appropriately mapped sensor data received aligned with the correct channel icon and/or text with each channel separated into sections on the display.
- 5.1.1.14. Error Message Display
  - 5.1.1.14.1. Error display if no Participant ID entered.
  - 5.1.1.14.2. Error display for access privileges denied for saving or accessing waveform data file(s).
  - 5.1.1.14.3. Error displayed on startup if configuration unavailable during startup.
- 5.1.2. File Management
  - 5.1.2.1. Recorded files are stored in default folder :  
C:\Users\Public\Documents\Recordings\ or other designated folder(s)
  - 5.1.2.2. Files are stored as CSV files.
  - 5.1.2.3. Files have the naming convention [Site][Task][Participant ID][Session] as selected and/or entered on the user interface.
- 5.1.3. Data Management / Conversion
  - 5.1.3.1. Invoke test message sent from USB to Bioamplifier. *[see Bioamplifier - Firmware Software Design for more information]*
    - 5.1.3.1.1. Automated Test message sent from the Visualizer System - Collect Software application upon initialization (Start Up) of the Visualizer system - Collect Software application.
    - 5.1.3.1.2. Manual Test message manually sent to trigger the Bioamplifier to respond with preconfigured EEG data in order to verify test results are received, the Bioamplifier is functioning properly and the Collect Software application is functioning properly.
  - 5.1.3.2. CSV Output files
    - 5.1.3.2.1. Recording output files [see section 5.2.2]. Recording file has 2 parts:Header and Body
    - 5.1.3.2.2. Header
      - JSON (JavaScript Object Notation) open-standard format.
      - 1 - Collect recording
        - Specification for recording file header
      - 2 - Channels
        - Channel description
        - 2.1 - Gain
          - ADC gain value
        - 2.2 - Index



- 0-based column index for this channels data in a row
- 2.3 - IsOn
  - Channel enabled
- 2.4 - label
  - Channel label
- 3 - Channel count
  - Number of channels in each row data
- 4 - Timestamp index
  - 0-based column index for timestamp data in a row
- 5 - Mark index
  - 0-based column index for mark data in a row
- 6 - Unit
  - Channel data unit (microvolt or ADC)
- 7 - Sample rate
  - Channel data sample rate in hz
- 8 - Collect version ( collect software application version)
  - Collect software version that generated this file
- 9 - Bioamplifier ID
  - Identifier for the EEG Bioamplifier hardware used
- 10 - Bioamplifier Version
  - EEG Bioamplifier Firmware SW Version
- 11 - Start record datetime
  - Date and time when the recording is started
- 12 - Site
  - Site that generated the recording
- 13 - Task
  - Task name performed during the recording
- 14 - Participant ID
  - Participant identifier for the recording
- 15 - Session
  - Session number for the recording

#### 5.1.3.2.3. Body

Formatted section separated from Header section using format

```
<CR><LF><CR-----><LF>
```

channelCount: [# of channels]

Channels:

```
{ gain: [#],  
  Index: [#],  
  isOn: [true]/[false],
```



```
Label: [string] e.g. "Fz",  
,  
{ etc... repeat for each channel....},  
{  
timestampIndex: [#],  
markIndex: [#],  
unit: [string],  
collectVersion: [#.#.#],  
bioampId: [string],  
bioampVersion: [string], }  
-----
```

Rows of comma-separated values.

Data in raw ADC output

Channel gain

One column of mark data

### 5.1.3.3. Configuration File presets

5.1.3.3.1. Configuration settings from first time use are stored in default Presets folder.

5.1.3.3.2. Saved Presets stored in executable root directory in [presets] folder (loaded automatically on next use).

5.1.3.3.2.1. Automatically updated via Collect Software application.

5.1.3.3.2.2. [preset name]\_preset.json

Channels:

```
{ gain: [#],  
  isBiasNOOn: [true]/[false],  
  isBiasPOOn: [true]/[false],  
  isOn: [true]/[false],  
  Label: [string] e.g. "Fz",  
},  
{ etc... repeat for each channel....},  
"Name": [name of preset] }
```

\* All possible channels listed with IsOn false for those not enabled.

### 5.1.4. Other

5.1.4.1. Keyboard and touchpad allow utilization of all software functionality.




- 5.1.4.2. WiFi disabled for utilization of Collect Software.
- 5.1.4.3. All software installed in root directory of executable.

## 5.1.5. Performance

- 5.1.5.1. Files are automatically saved during the selection of the Stop Recording command button option.
- 5.1.5.2. EEG Data is continuously received from the connected USB port and captured automatically during the 'Recording' process.
- 5.1.5.3. The Send Test Command Button option is available for sending test messages to the Bioamplifier in order to verify connectivity and capture of configured test channel results.

## 5.2. Bioamplifier

Fig 1. Bioamplifier PCB



### 5.2.1. Firmware Software Design

Internal and external drivers utilized for managing firmware response.

#### 5.2.1.1. EEG Data

- 5.2.1.1.1. EEG Data is continuously streamed on the DATA port in ASCII format.
- 5.2.1.1.2. Stimulus Mark data is streamed on the STIM port in ASCII format.
- 5.2.1.1.3. Channels comma separated and encoded as signed 24bit value for each. - handling up to 32 channels of input data.
- 5.2.1.1.4. Implement command line interpreter (CLI). command to CLI are terminated with carriage return.
- 5.2.1.1.5. EEG Data collection, Blinking LED's, and serial communication are managed through a main loop.
- 5.2.1.1.6. EEG data format:



SEQUENCE NUMBER	CH 1	CH 2	...	CH32	MARK	CR/LF
--------------------	------	------	-----	------	------	-------

- 5.2.1.1.6.1. EEG Data prefixed by 'DATA'
- 5.2.1.1.6.2. Sequence number: Auto incrementing 32bit value
- 5.2.1.1.6.3. Mark: user -selected data value
- 5.2.1.1.6.4. CR/LF terminates each line.
- 5.2.1.1.7. Event Driven or Sent Commands
  - 5.2.1.1.7.1. Echo STIM Port commands
  - 5.2.1.1.7.2. Multiplex DATA Port commands with EEG Data and prefixed with 'CLI'.
  - 5.2.1.1.7.3. Additional Commands
    - EEG Data ready
    - Serial Comm
    - Power on / off
- 5.2.1.2. Test Command
  - 5.2.1.2.1. Data message sent from connected Visualizer System - Collect Software Application through serial USB port on Laptop PC, connected USB communication cord, to the DATA USB port on the Bioamplifier.
  - 5.2.1.2.2. Register settings are configured to generate internal amplitude and frequency settings and output to the calling source (the Visualizer System - Collect Software Application)
    - Input - Test signal to trigger Bioamplifier response
    - Output - Test response with internally generated Amplitude (uV) and Frequency (Hz)
- 5.2.2. Hardware Design
  - 5.2.2.1. Classification
    - 5.2.2.1.1. FDA Class II Medical Device
  - 5.2.2.2. Overall Layout
    - 5.2.2.2.1. Enclosure - 155mm (W) x 148mm (D) x 57mm (H), 3.3mm thickness
    - 5.2.2.2.2. Unit Weight < 1 kg
  - 5.2.2.3. Operation Mode
    - 5.2.2.3.1. Continuous Operation Mode
  - 5.2.2.4. Temperature Limits
    - 5.2.2.4.1. Operating Temperature 0 to 30°C

For Internal Use Only. This document contains confidential, proprietary information.  
It may not be copied or reproduced without permission.



- 5.2.2.4.2. Non-Operating Temperature -20 to 60°C
- 5.2.2.5. Data Acquisition (DAQ)
  - <http://www.ti.com/lit/ds/symlink/ads1299.pdf>
  - 5.2.2.5.1. (4) IC ADS1299 - 8 channel bioamplifier (Texas Instruments)
  - 5.2.2.5.2. (1) ADS1299 serves as the master (U1A) and provides the conversion clock to the additional (3) units for synchronizing as well as master bias generator.
  - 5.2.2.5.3. Parallel data configuration (integration).
  - 5.2.2.5.4. Ground bias generator.
  - 5.2.2.5.5. Built in lead-off detection.
  - 5.2.2.5.6. 24 bit conversion
  - 5.2.2.5.7. Channel Inputs
    - 32 Total EEG channels
    - Single-ended inputs for channels 1-24
    - Single-ended or differential inputs for channels 25-32
  - 5.2.2.5.8. 250 samples / sec Sample Rate
  - 5.2.2.5.9. MCU communication
    - 4 wire SPI @ clock speed of 4MHz
    - ADS1299 sends signals that lower for conversion readiness from the chips.
    - Master IC ADS1299 signal connection with MCU
    - Chip selects are asserted to read the data from each of the ADS1299s in turn.
  - 5.2.2.5.10. Consolidates Input EEG data with received STIM port data.
- 5.2.2.6. Microcontroller (MCU)
  - 5.2.2.6.1. 64 pin QFP package Kinetis (NXP/Freescale) KL27Z256
  - 5.2.2.6.2. Arm Cortex M0 processor running at 48 MHz with 256K flash memory.
- 5.2.2.7. Input Pins
  - (4) input pins - hardware revision (via pull-down resistors)
  - (4) pins reserved
- 5.2.2.8. Power mode
  - Full-power mode (default)
  - Low-power mode (optional)
- 5.2.2.9. Test and Programming
  - Test via dedicated JTAG port implemented via TAG connect cable.
  - Programming via DATA USB port using built-in bootloader and proprietary tool.



- 5.2.2.10. Power Supply (Detachable)
  - 5.2.2.10.1. Via (2) USB connectors (either may power board)
  - 5.2.2.10.2. Allow powering from external source for 90 minutes minimum length.
  - 5.2.2.10.3. USB Voltage - diode-ORed
    - Powers isolated DC to DC Converter (PS1)
    - 4kv isolation and patient leakage current protection
  - 5.2.2.10.4. DC to DC Converter
    - 5V output linear regulated down to 3.3V (low noise floor)
  - 5.2.2.10.5. 3.3V down voltage:
    - Powers MCU and digital logic.
    - Powers charge pump based converter that provides the +/-2.5v rails to the DAQ ADS1299(s).
  - 5.2.2.10.6. 2.5V reference power
    - Powers MCU Analog section
    - Provide precise reference for measurements
  - 5.2.2.10.7. Power supply Source
    - Input - 100-240 Vac 50/60 Hz 1.8A (Dell Adapter HAI30PM130)
    - Output - 19.5Vdc, 7.67A
  - 5.2.2.10.8. Overvoltage Category - OVC II
- 5.2.2.11. Coin cell battery on PCB
  - 5.2.2.11.1. CR2332 - Lithium Magnesium , 3V, 225mAH, 60°C
- 5.2.2.12. USB Ports
  - 5.2.2.12.1. (2) USB Ports for PCB Connection to (2) PC's.
    - Port 1 - DATA Port - Transfer streaming EEG data collected.
    - Port 2 - STIM Port - Board control and sets marks for the corresponding streaming EEG data collected.
  - 5.2.2.12.2. PC Side - FTDI232 USB to serial IC
  - 5.2.2.12.3. MCU Side - optical isolators for patient protection.
  - 5.2.2.12.4. Rated at 4kv - max data rate 10mbps
  - 5.2.2.12.5. USB Type B, 4 contacts, 1.5A, UL94V-0, Plastic housing, 105°C
- 5.2.2.13. Real Time Clock
  - 5.2.2.13.1. Communicates over I2C with the MCU
  - 5.2.2.13.2. Battery backup of real time clock.
- 5.2.2.14. Micro SD
  - 5.2.2.14.1. Connected to a second SPI port on the MCU
- 5.2.2.15. PCB





- 5.2.2.15.1. Circuitboard - 155mm x 75mm, 1.6mm thickness, ULV-0, 130°C
- 5.2.2.15.2. 6 layer PCB board  
Layer Stack signal > ground plane > signal > signal > power connections > signal
- 5.2.2.15.3. Minimum line/space is 6/6 and minimum via size is 6/16.
- 5.2.2.15.4. Board material is plain FR4 and finish is ENIG.
- 5.2.2.15.5. Board is separated into (3) sections (PC-non isolated, Digital, Analog)  
PC-non isolated - USB Ports, half of the optical isolators, and the DC to DC converters (Creepage and clearance - 8mm)  
Digital - Power Supplies, MCU and expansion header  
Analog - Protects the uV EEG Signals from digital noise (split ground plane). An RF shield for encapsulating the ADS1299.

Fig 2 Bioamplifier PCB (with sections)





Fig 3. Bioamplifier PCB (with split ground plane)



## 5.2.2.16. Unit connector

- 5.2.2.16.1. PCB D-Sub Connector, receptacle type, 68 pos, 1.27mm pitch, 1A max contact current, UL 94V-0, 105°C
- 5.2.2.16.2. Allows connectivity to Adapter (Headset) attachment - male to female connector setting.

## 5.2.2.17. Gasket

- 5.2.2.17.1. O-ring, vulcanized EPDM rubber
- 5.2.2.17.2. 19mm ring diameter, 1.78mm thickness, 150°C

## 5.2.2.18. LED's (Light Emitting Diodes)

- 5.2.2.18.1. LED for display 'On' during powering of device.
- 5.2.2.18.2. LED for port connection (2) STIM and DATA USB Ports
- 5.2.2.18.3. LED (4) for activation, programming / test and data communication.

## 5.2.2.19. Parts Reference

Part Number	Part Type	Description
710-03064-01	710-EM Assembly	ASSY, M20 BIOAMP, BUTTER DISH, AMBER
690-13322-01	690-Mech Part	MECH, BUTTER AMP TOP, AMBER
690-13321-01	690-Mech Part	MECH, BUTTER AMP BOTTOM
690-13335-01	690-Mech Part	MECH, LIGHTPIPE, BUTTER
650-02461-01	650-PCBA	PCBA, LUCHADOR, REV 3.2b, AMBER



690-13336-01	690-Mech Part	MECH, GASKET, LIGHTPIPE, BUTTER, AMBER
690-13337-01	690-Mech Part	MECH, BOOT, CONNECTOR, BUTTER
250-06431-00	250-Hardware Part	HW, SCREW, THREAD FORMING, TORX RDH, M2.63 x 1.06 THD x 6 MM LG, 18-8 SST, AMBER
250-06430-00	250-Hardware Part	HW, SCREW, THREAD FORMING, TORX RDH, M3.63 x 1.34 THD x 10 MM LG, 18-8 SST, AMBER
690-13338-01	690-Mech Part	MECH, FOOT, BUTTER
690-13950-01	690-Mech Part	MECH, O-RING, 191.26MM ID X 1.78THK, AMBER
690-13430-01	690-Mech Part	MECH, GASKET, USB, BUTTER, AMBER
690-13885-01	690-Mech Part	MECH, TAPE, VHB RING, AMBER
961-00019-01	952-Label	LABEL, INVESTIGATION, REGULATORY, AMBER

### 5.2.3. Performance

- 5.2.3.1. MCU must be fast with high availability and simple to develop with product 'tool chain'.
- 5.2.3.2. IC (LM27762) - built-in linear regulator to lessen noise output.
- 5.2.3.3. Proprietary firmware drivers to increase performance and lower latency.

### 5.2.4. Safety


- 5.2.4.1. Excessive current and voltage protected by circuit design.

### 5.2.5. Label

- 5.2.5.1. Philo M20 Bioamp device shall be labeled per GOO-351-CSD-961-00019-01-DWG .
- 5.2.5.2. Serial Number format is as follows: BYYMMDD##, where YY are the last two digits of the current year, MM is the current month, DD is the current day, and ## are the number in sequences of devices made on that day (starting with 01).



Fig. 4 Bioamplifier Label



## 5.3. Adapter (Headset Unit) w/ Sensor Strip

### 5.3.1. Hardware Design

#### 5.3.1.1. Classification

##### 5.3.1.1.1. FDA Class II Medical Device

#### 5.3.1.2. Overall Layout (Adapter)

##### 5.3.1.2.1. Enclosure - 137mm (L) x 134mm (W) x 49mm (H), 3.3mm thickness

##### 5.3.1.2.2. Unit Weight < 1 kg

#### 5.3.1.3. Operation Mode

##### 5.3.1.3.1. Continuous Operation Mode

#### 5.3.1.4. Temperature Limits

Operating Temperature 0 to 30°C

##### 5.3.1.4.1. Non-Operating Temperature -20 to 60

#### 5.3.1.5. PCB


##### 5.3.1.5.1. Circuitboard - 94mm x 75mm, 1.6mm thickness, ULV-0, 130°C

##### 5.3.1.5.2. 100uA current source for each EEG channel and powers the electrode (operational electrode in bootstrap mode)



- 5.3.1.5.3. uV signals from operational amplifier are high-pass filtered at .1Hz buffered and low-pass filtered at 40Hz with both signals passed to Bioamplifier.
- 5.3.1.5.4. Reference electrode processed as EEG electrode.
- 5.3.1.5.5. Bias electrode is pass-through without processing.
- 5.3.1.5.6. PCB Board - 4 layer board

Fig. 5 Headset Adaptor PCB



- 5.3.1.6. Unit Connector
  - 5.3.1.6.1. PCB D-Sub Connector, plug type, 68 pos, 1.27mm thickness, 1A max contact current, UL 94V-0, 150°C (SCSCI connector)
  - 5.3.1.6.2. Allows secure connectivity to Bioamplifier unit - male to female connector setting.
- 5.3.1.7. Sensor Strip - Strip Body
  - 5.3.1.7.1. Synthetic leather - 225mm (L) x 37mm (W), 2.6mm thickness
- 5.3.1.8. Sensors
  - 5.3.1.8.1. (3) or (4) - 19mm diameter metal base, 1.9mm in diameter 6.6mm high metal, 8 equal distances distributes around base parameter.
  - 5.3.1.8.2. Bias and Reference ear clip attachments.
  - 5.3.1.8.3. Transfer signals collected to the adapter.
  - 5.3.1.8.4. High impedance sensors.
- 5.3.1.9. Headset - Cap








- 5.3.1.9.1. The # of sensing electrodes (sensors) are placed in the Fz, Cz, and Pz positions (Oz optional), Two additional electrodes are clipped to the earlobes at A1 and A2.
  - 5.3.1.9.1.1. Caps has holes for attaching 3 sensor or 4 sensor headset attachments.
- 5.3.1.9.2. Chinstrap attachment
- 5.3.1.9.3. Flexible material with industry compliant holes positioned for sensors attachment. Flexibility includes ability to apply the Cap with Sensor strip within 3 minutes on participant and will fit a participant with a head diameter of 52-64cm. Verify ability to remove the headset and sensors from the participant in less than 1 minute.
- 5.3.1.9.4. Durable to enviroid cleaning agent applied up to 9 times to Cap and remain functioning for use. Durability includes ability to use at least 10 times for participant sessions. The sensors and associated wires are protected by the headset.
- 5.3.1.9.5. The sensors protrude up to 20mm from the headset to allow for hair thickness.
- 5.3.1.9.6. Headset / Adaptor cable length - 50-250cm.
- 5.3.1.9.7. The Textile is 'light weight' and may be worn under reasonable conditions and allow for signal collection for 60 minutes minimal.
- 5.3.1.9.8. 2 clip attachments for ear (Bias and Reference sensors)
- 5.3.1.10. Power Supply
  - 5.3.1.10.1. Allow powering from external source for 90 minutes minimum length.
- 5.3.1.11. Communication
  - 5.3.1.11.1. EEG Signals collected from Sensors are transmitted to the Adapter.
  - 5.3.1.11.2. Active Sensor Signals from Adapter are then transferred to the Bioamplifier.
- 5.3.1.12. Cleaning
  - 5.3.1.12.1. Allows wiping down with a soft dry cloth or germicidal disposable wipe.
- 5.3.1.13. Label
  - 5.3.1.13.1. Philo M20 Headset - Adapter device shall be labeled per PL-0322-00
  - 5.3.1.13.2. Serial Number format is as follows: BYYMMDD##, where YY are the last two digits of the current year, MM is the



# EEG Product Design Specification

xxx-xxxxx-00  
Rev 1

current month, DD is the current day, and ## are the number in sequences of devices made on that day (starting with 01).

		<b>CAUTION--Investigational device. Limited by Federal (or United States) law to investigational use.</b> <b>FCC Statement: This device has not been authorized as required by the rules of the Federal Communications Commission. This device is not, and may not be, offered for sale or lease, or sold or leased, until authorization is obtained.</b>
<b>SN</b>	H19031104	
<b>REF</b>	BOM 0210-01	
<b>QTY</b>	1 6-Sensor Headset	
	Google LLC, Mountain View, California 94043	
	Read instructions before use	
	Type BF applied part	
Contact: M20System@google.com		PL-0322-00

## 5.3.1.14. Parts Reference

Part Number	Part Type	Description
690-13313-01		MECH BUTTER ADAPTER BOTTOM
690-13320-01		MECH BUTTER ADAPTER TOP
250-06431-00		HW SCREW, THREAD FORMING,TORX RDH, M2.63 x 1.06 THD x 6 MM LG, 18-8 SST, AMBER
650-02489-01		PCBA, A.E BREAKOUT, AMBER
690-13337-01		MECH, BOOT, CONNECTOR, BUTTER
250-06430-00		HW, SCREW, THREAD FORMING, TORX RDH, M3.63 x 1.34 THD x 10 MM LG, 18-8 SST, AMBER
690-13334-01		MECH, GROMMET, BUTTER, AMBER
690-14175-01		MECH, O-RING, 149.35MM ID, 1.78 THK, AMBER
710-03081-01		ASSY, ELEC-MECH, M20 ADAPTER, BUTTERSTICK, AMBER
725001034-01		ASSY, HARNESS, WIRE HARNESS PLUS SENSOR STRIP, M20, AMBER
210-04100-00		CONN, CLIP, G.SAHARA, ACTIVE ELECTRODE, AMBER

For Internal Use Only. This document contains confidential, proprietary information.  
It may not be copied or reproduced without permission.





210-04098-00		CONN, REF CLIP, G.ACTIVE ELECTRODE SYSTEM, AMBER
690-14183-01		MECH, SENSOR LEATHER STRIP, M20, BOTTOM, GRAY, AMBER
210-04101-00		CONN, GROUND CLIP, G.ACTIVE ELECTRODE SYSTEM, AMBER
690-14184-01		MECH, SENSOR LEATHER STRIP, M20, TOP, YELLOW, AMBER
250-06645-00		HW, TAPE, CLEAR, ADHESIVE DOUBLE SIDED, 3/8IN WD, AMBER
250-06657-00		HW, CABLE TIE, 4 IN LG, 10 LB BREAK STRENGTH, NYLON, OFF-WHITE, AMBER
250-06641-00		HW, SLEEVING, FLEXO, 3/8 IN, 3/16 IN to 1/2 IN EXP, 500FT, WHT, AMBER
PN-0057-00	SAHARA clip	SAHARA active electrode clip connector
PN-0058-00	ACTIVEclipREF	Reference clip connector
PN-0060-00	ACTIVEclipGND	Ground clip connector
2535-03	Adhesive Tape	Tanner's Bond Adhesive Tape
7120K12	Cable Tie	Narrow, 4" Long, 10 lbs. Breaking Strength, Off-White
PN-0101-00	Flexo PET	Flexo PET (PTN)

## 5.4. USB Cables

5.4.1. USB 2.0 A (male)/ B(male), 28A WG data wires, 24A WG power wires, 1.8m length.

5.4.1.1. Cable for connectivity from Bioamplifier to Visualizer System.

5.4.1.2. Cable for connectivity from Stimulus System to Bioamplifier.

## 5.5. Laptop Computer

5.5.1. Execution of the Visualizer System (Collect software application)

5.5.1.1. Icon available on the desktop for selection

5.5.1.2. Directory for saving recordings - desktop shortcut

5.5.1.3. 3 Operating System Logins

5.5.1.3.1. Lab Technician

For Internal Use Only. This document contains confidential, proprietary information.  
It may not be copied or reproduced without permission.





Utilized by site staff for execution of M20 site tests

5.5.1.3.1.1. Limited operating system access

5.5.1.3.1.2. Full access to applications

5.5.1.3.2. Rocket Turtle

Utilized by X-Team staff for managing of M20 site studies.

5.5.1.3.2.1. Increased operating system access

5.5.1.3.2.2. Full access to M20 and most system applications

5.5.1.3.3. Administrative

Restricted to X-Team staff administrators

5.5.1.3.3.1. Full operating system access

5.5.1.3.3.2. Full access to all system application

5.5.1.4. Collect Software application has default user permissions for read / write / delete / execute within directory as well as to output file directory for saved recordings.

5.5.1.5. Collect Software application has full Port Send / Receive access.

5.5.1.6. Collect Software application maximizes application screen resolution to full Laptop screen resolution.

5.5.1.7. System displays the current date and time on the bottom right hand side of the screen. HH:MM AM/PM & MM/DD/YYYY

## 5.5.2. Software Specification

5.5.2.1. Operating system software - Microsoft Windows 10 Professional 64 bit

5.5.2.2. Collect software application ver. 2.2.3

5.5.2.3. Disk encryption - BitLocker full disk encryption

5.5.2.4. Google Chrome ver. 75+ 64 bit

5.5.2.5. FTDI Driver USB Serial driver supported on Windows 10 Professional

5.5.2.6. DirectX ver 9+

5.5.2.7. Installation Dependencies

5.5.2.7.1. Bower 1.8.8

5.5.2.7.2. Debug 4.1.1

5.5.2.7.3. Electron-Settings 3.2.0

5.5.2.7.4. Express 4.16.2

5.5.2.7.5. Gasp 2.0.2

5.5.2.7.6. Opn 5.4.0

5.5.2.7.7. Serial Port 6.2.1

5.5.2.7.8. Sudo Prompt 8.2.3

5.5.2.7.9. WinReg 1.2.4

5.5.2.7.10. WS 6.0.0

5.5.2.8. Default Presets

For Internal Use Only. This document contains confidential, proprietary information.  
It may not be copied or reproduced without permission.



5.5.2.8.1. [executable install root directory] > Resources > App > presets

## 5.5.3. Hardware Specification

5.5.3.1. Processor - i7+, 2.6GHz+ minimum

5.5.3.2. Harddrive - 500GB+ HD

5.5.3.3. Memory - 32GB+ RAM Memory

5.5.3.4. Battery

5.5.3.4.1. 56 WHr+

5.5.3.4.2. 19Vdc input

5.5.3.4.3. Battery will allow 90 minutes+ continuous power without recharge.

5.5.3.5. Keyboard for data entry and selection/navigation.

5.5.3.6. Monitor - 15"+ (1x15x10)

5.5.3.7. Graphics - 59Hz, Intel UHD Graphics 630

5.5.3.8. USB Ports

5.5.3.8.1. Bioamplifier - Serial over USB 2.0 Port

5.5.3.8.2. External Mouse port - USB 2.0 Type A port

5.5.3.9. Default Inactive Lock setting set for 60 minutes. 5.5.3.9